



Technical Resources

Protocols

Vectors

Tools

Citations

FAQs

eNotes

Promega Notes

Cell Notes

Guides

MSDS

Publications

Contact Us

Literature Request

New Products

All Applications

FAQs

RNasin® Plus RNase Inhibitor

See Technical Resources for more information.

1. What is RNasin® Plus RNase Inhibitor?
2. Will RNasin® Plus RNase Inhibitor work better than native or Recombinant RNase Ribonuclease Inhibitor?
3. How is RNasin® Plus RNase Inhibitor more resistant to oxidative stress than recombinant native RNasin®?
4. What are the characteristics of RNasin® Plus RNase Inhibitor?
5. How can I use the characteristics of RNasin® Plus RNase Inhibitor to better protect my RNA template?
6. How do the characteristics of RNasin® Plus RNase Inhibitor relate to higher temperature reverse transcriptase (RT) reactions?
7. What applications are compatible with RNasin® Plus RNase Inhibitor?

1. What is RNasin® Plus RNase Inhibitor?


RNasin® Plus RNase Inhibitor is a recombinant mammalian RNase inhibitor expressed as a soluble protein in *E. coli*. Through natural amino acid diversity RNasin® Plus RNase Inhibitor has increased resistance to oxidation when compared to the human protein.

2. Will RNasin® Plus RNase Inhibitor work better than native or Recombinant RNase Ribonuclease Inhibitor?

RNasin® Plus RNase Inhibitor is naturally more resistant to oxidative stress, increasing overall efficacy. As an added feature RNasin® Plus displays continued inhibition at higher temperatures. Also, RNasin® Plus is expressed by *E. coli* as a soluble protein, allowing easy purification by a combination of ion exchange and hydrophobic interaction chromatography. No direct affinity chromatography is required. This new process yields >90% pure protein with no *E. coli* RNase carryover. However, the mechanism of action remains the same for RNasin® Plus and native or Recombinant RNasin® Inhibitor. The mechanism is an inhibition of eukaryotic RNases via stoichiometric 1:1 noncovalent binding of the RNasin® Plus RNase Inhibitor to an RNase.

3. How is RNasin® Plus RNase Inhibitor more resistant to oxidative stress than native or recombinant RNasin®?

Two cysteines in the human protein have been identified as especially sensitive to oxidation and react by forming a disulfide that can block the active site of the inhibitor (1). RNasin® Plus RNase Inhibitor, through natural amino acid diversity, lacks the ability to form a blocking disulfide.



4. What are the characteristics of RNasin® Plus RNase Inhibitor?

During development of the RNasin® Plus RNase Inhibitor Promega scientists discovered continued inhibition of RNases even above the normal denaturation temperature of the RNasin® Plus molecule. A mixture of RNasin® Plus and a pure RNase, like RNase A heated to at least 70°C for 15 minutes, and the RNase A activity does not return upon cooling to normal temperatures, such as those used in the RT step of RT-PCR. It has been demonstrated to work in this manner with a complex mixture of RNases in a rat liver protein extract (Sigma Cat.# L1380). Rat liver is known to contain ribonucleases (2). No detectable RNase activity, as determined by RT-PCR, is observed when a mixture of rat liver RNases and RNasin® Plus RNase Inhibitor is heated for 15 minutes followed by the addition of 100ng or 10ng of template RNA and incubation for an additional hour at 37°C.

5. How can I use the characteristics of RNasin® Plus RNase Inhibitor to protect my RNA template?

Many protocols, including those for ImProm-II™ Reverse Transcription System require an initial thermal denaturation of the RNA template of interest in the presence of reverse transcription primers for 5–10 minutes at 70°C followed by a quick chill step denatures secondary structure in the RNA template, allowing greater sensitivity for PCR. In light of the new activities identified for RNasin® Plus RNase Inhibitor, it can now be added at this step to protect the RNA template during thermal denaturation. RNases that were present during the thermal denaturation will be inactivated; hence, no more RNase inhibitor should be added during full RT reaction assembly in the event additional exogenous RNases are inadvertently added to the reaction from pipette components or other sources.

6. How do the characteristics of RNasin® Plus RNase Inhibitor relate to high temperature reverse transcriptase (RT) reactions?

The characteristics of RNasin® Plus RNase Inhibitor allow you to set up your high temperature first-strand synthesis reactions and take them to reverse transcriptase reaction temperatures above 50°C. This gives researchers RNase protection when transcribing RNA templates with high secondary structure.

7. What applications are compatible with RNasin® Plus RNase Inhibitor?

RNasin® Plus RNase Inhibitor has been tested in RT-PCR and is compatible with such as AMV, M-MLV and ImProm-II™ Reverse Transcriptases or Taq and T7 DNA Polymerases. RNasin® Plus RNase Inhibitor has also been tested and found to be compatible with quantitative, real-time RT-PCR reactions in a TaqMan® Assay. The new inhibitor is also compatible with the Riboprobe® System for in vitro transcription using the T3, T7 RNA Polymerase. RNasin® Plus RNase Inhibitor can also be used with Wheat Germ and Rabbit Reticulocyte Lysate for in vitro translation from an RNA template, as in the TNT® Wheat Germ and TNT® Reticulocyte Lysate System for coupled in vitro transcription/translation.

References

1. Kim, B.M., Schultz, L.W. and Raines, R.T. (1999) Variants of ribonuclease inhibitor resist oxidation. *Protein Sci.* **8**, 430–4.
2. Zhao, W. et al. (1998) Ribonucleases from rat and bovine liver: Purification, spectroscopic and structural characterization. *Biochim. Biophys. Acta* **1384**, 55–65.